

Micro Screw System & Micro Screw System Basic

,Olsberger Konzept' by Prof. Dr. Fouad Khoury

The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.



BMSBA
Micro Screw System Basic



BMS00
Micro Screw System

At a glance

- Safe fixation and stabilization of cortical bone grafts
- Osteosynthesis screws with diameters of 1.0 and 1.2 mm
- Osteosynthesis screws made of surgical stainless steel for high stability
- Self-tapping thread flanks allow the osteosynthesis screws to be screwed in without creating a gliding hole

Our sales force would be pleased to answer your questions: +49 2131 2012-303

Meisinger since 1888

Micro Screw System Kits

Fax for ordering +49 2131 2012-222

Figure	Shank	Size	Min.-/External-Ø/ Length [mm]	Piece / Kit	Packing unit for refill	Amount	
	BMSBA	-	-	-	1	1	
	MSPB1	204	008	- / 0.8 / 14.0	1	5	
	MSPB2	204	010	- / 1.0 / 14.0	1	5	
	202RF	206	008	- / 0.8 / 14.0	1	5	
	202RF	206	010	- / 1.0 / 14.0	1	5	
	MSS10 060	-	-	0.6 / 1.0 / 6.0	3	3	
	MSS10 080	-	-	0.6 / 1.0 / 8.0	3	3	
	MSS10 100	-	-	0.6 / 1.0 / 10.0	3	3	
	MSS12 060	-	-	0.7 / 1.2 / 6.0	3	3	
	MSS12 080	-	-	0.7 / 1.2 / 8.0	3	3	
	MSS12 100	-	-	0.7 / 1.2 / 10.0	3	3	
	MSSDE	-	-	-	1	1	
	MSSDW	204	-	-	1	1	
	BMS00	-	-	-	1	1	
	MSPB1	204	008	- / 0.8 / 14.0	1	5	
	MSPB2	204	010	- / 1.0 / 14.0	1	5	
	202RF	206	008	- / 0.8 / 14.0	1	5	
	202RF	206	010	- / 1.0 / 14.0	1	5	
	MSS10 040	-	-	0.6 / 1.0 / 4.0	2	3	
	MSS10 060	-	-	0.6 / 1.0 / 6.0	3	3	
	MSS10 080	-	-	0.6 / 1.0 / 8.0	3	3	
	MSS10 100	-	-	0.6 / 1.0 / 10.0	3	3	
	MSS10 120	-	-	0.6 / 1.0 / 12.0	2	3	
	MSS10 140	-	-	0.6 / 1.0 / 14.0	2	3	
	MSS12 040	-	-	0.7 / 1.2 / 4.0	1	3	
	MSS12 060	-	-	0.7 / 1.2 / 6.0	2	3	
	MSS12 080	-	-	0.7 / 1.2 / 8.0	2	3	
	MSS12 100	-	-	0.7 / 1.2 / 10.0	2	3	
	MSS12 120	-	-	0.7 / 1.2 / 12.0	1	3	
	MSS12 140	-	-	0.7 / 1.2 / 14.0	1	3	
	MSSDM	-	-	-	1	1	
	MSGSD	-	-	-	1	1	
	MSSDW	204	-	-	1	1	

Hager & Meisinger GmbH

Hansestr. 10 • 41468 Neuss • Germany
Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222
info@meisinger.de • www.meisinger.de

Meisinger USA, L.L.C.

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
Toll free: +1 (866) 634-7464
info@meisingerusa.com • www.meisingerusa.com

Your delivery address:

